

FFP2 KN95 Four-Layer Filtration Protective Mask

SPECIFICATION

PT-103A



Multiple
protection



KN95
grade



Light &
breathable



3D fitting



Skin-friendly



Bacterial
filtration




Up to 95% filtration efficiency
KN95 Standard


4-layer
filtration protection

Features:

- 4 layer lining
- Up to 95% filtration efficiency
- Air pollution resistance (PM 2.5)
- Comfortable 3D flex design
- Reusable & eco-friendly
- Suitable for adults and children



Skin-friendly, soft non-woven fabric
low sensitivity and no irritation



Meltblown Fabric, key filtration material

KN95 filtration efficiency reaches 95%,
effectively protecting respiratory health



Adjustable elastic ear hanging rope




Respiratory resistance	KN95 filter effect	Formaldehyde content
< 80 pa	95 %	0 %

Specifications

Product Name	FFP2/FFP3 KN95 Face Mask
Aerosol Type	Non-Oil
Assigned Protection Factor	10pcs
Number per Box	50pcs
Boxes per Case	20box
Case Quantity	1000pcs
Class	Comfort
Exhalation Valve	No
Features	Advanced Electrostatic Media
Product Color	White
Product Type	Dust and other Particles
Recommended Industry	Commercial Buildings,Construction,Design & Construction, Food Processing,Food Safety,General Manufacturing, Heavy Infrastructure,Mining,Oil & Gsa,Transportation
Standards/Approvals	KN95/ CE EN 149:2001+A1
Expiration Date	2 years

Product Dimensions



Packaging Details

Product



1 pc

Protective Bag



5pcs/bag

Box



10 bags
50 pcs/box

Case




20 boxes
1000 pcs/case

Packaging Dimensions & Weight

Box Dimension and Weight

Size: 5.12"L x 4.92"W x 7.28"H

Weight: 0.83lbs



Case Dimension and Weight

Size: 25.39"L x 10.04"W x 0.29"H

Weight: 17.64lbs



User Guide



Step 1.

With the outer nose-clip facing up, place the mask over the face to cover both nose and mouth.

Step 2.

Pull the left and right ear straps around each ear for a comfortable fit.

Step 3.

Using both hands, press the outer nose-clip inward to comfortably form to and support the nose bridge.

Step 4.

Check for air leakage around the face by breathing in-and-out and continue to adjust the mask to achieve a proper fit.



Fiscal Year 2020

CERTIFICATE OF FDA REGISTRATION

This certifies that:

SHENZHEN PEARLTECK CO.,LTD.

**Room 0906,Section B,Rongchao Binhai Building,South of Xinghua Road,Xin'an
Sub-district,Bao'an District, Shenzhen, Guangdong, 518101, CHINA**

has completed the FDA Establishment Registration and Device Listing with the US Food & Drug Administration, through UCL-REG SERVICE INC.

Owner/Operator Number: 10065655

Listing No.
D381936

Product Code:
QKR

Device Name:

Face mask (except N95 respirator) for general
public/healthcare personnel per IIE guidance
Face Mask PT-103A, PT-103X "X" stands for numbers
from 0 to 9, letters from A to Z

UCL-REG SERVICE INC. will confirm that such registration remains effective upon request and presentation of this certificate until the end of the calendar year stated above, unless said registration is terminated after issuance of this certificate. UCL-REG SERVICE INC. makes no other representations or warranties, nor does this certificate make any representations or warranties to any person or entity other than the named certificate holder, for whose sole benefit it is issued. This certificate does not denote endorsement or approval of the certificate-holder's device or establishment by the U.S. Food and Drug Administration. UCL-REG SERVICE INC. assumes no liability to any person or entity in connection with the foregoing.

Pursuant to 21 CFR 807.39, "Registration of a device establishment or assignment of a registration number does not in any way denote approval of the establishment or its products. Any representation that creates an impression of official approval because of registration or possession of a registration number is misleading and constitutes misbranding."

The U.S. Food and Drug Administration does not issue a certificate of registration, nor does the U.S. Food and Drug Administration recognize a certificate of registration, UCL-REG SERVICE INC. is not affiliated with the U.S. Food and Drug Administration.



UCL-REGSERVICE INC.

602 ROCKWOOD ROAD,WILMINGTON,

NEW CASTLE DE 19802 USA

UCL-REG SERVICE INC.

Cert. No.: M20738-2

Issued Date: 30 March 2020

Expiration Date: 31 December 2020

Certificate of Compliance



No. 0D200405R.SPC0S07

Technical Construction File no. SPC-20000402

Certificate's Holder: Shenzhen Pearlteck Co., Ltd
Room 0906, Section B, Rongchao Binhai Building,
South of Xinghua Road, Xin'an Sub-district,
Shenzhen, Guangdong, 518101, China)

Manufacturer: Dongguan Meikang Environment
Intelligent Technology Co. Ltd
No.4, Shaliba Street, Chang'an Town,
Dongguan City, Guangdong Province, China

**Certification ECM
Mark:**



Product: KN95 Face Mask
Model(s): PT-103X ('X' stands for numbers from 0 to
9, letters from A to Z)
Verification to: Standard:
EN 149:2001+A1:2009

related to CE Directive(s):
R 2016/425 (Personal Protective Equipment)

Remark: This document has been issued on a voluntary basis and upon request of the manufacturer. It is our opinion that the technical documentation received from the manufacturer is satisfactory for the requirements of the ECM Certification Mark. The conformity mark above can be affixed on the products accordingly to the ECM regulation about its release and its use.

Additional information and clarification about the Marking:



The manufacturer is responsible for the CE Marking process, and if necessary, must refer to a Notified Body. This document has been issued on the basis of the regulation on ECM Voluntary Mark for the certification of products. RG01_ECM rev.3 available at: www.entecerma.it

Issuance date: 05 April 2020

Expiry date: 04 April 2025

Reviewer
Technical expert
Amanda Payne



Approver
ECM Service Director
Luca Bedonni



Ente Certificazione Macchine Srl

Via Ca' Bella, 243 – Loc. Castello di Serravalle – 40053 Valsamoggia (BO) - ITALY
☎ +39 051 6705141 📠 +39 051 6705156 ✉ info@entecerma.it 🌐 www.entecerma.it